

REF AL0604



IVD



Σ 100

Version: 03, 13/06/2022



Axon Biotech GmbH
Im Oberfeld 2, DE-94491 Hengersberg
mail@medco.eu, www.medco.eu
A Company of Axonlab

Distributor:

Axon Lab AG
www.axonlab.com

NaCl 0.85%-solution for Kirby-Bauer Method

Intended use

The NaCl 0.85% solution for Kirby-Bauer is a ready-to-use solution for the agar diffusion test. It is an in vitro diagnostics product.

The ready-to-use test tubes containing 2 ml of the 0.85% NaCl solution can be used directly to prepare the inoculum of the desired test organism for the agar diffusion method.

Contents

100 Tubes NaCl 0.85% à 2 ml

Composition in %	
Sodium Chloride	0.85

Warning

The ready-to-use NaCl solution is intended for professional, single use only. Do not use if evidence of microbial contamination, discoloration, ruptures, or other signs of deterioration; ensure proper product disposal.

Storage

The NaCl 0.85% solution for Kirby-Bauer method is stored at 15°C to 30°C. Do not freeze. Do not use after the expiration date specified on the package. After opening the shipping package, unopened tubes can be used up until the expiry date if storage conditions have been met.

Quality control

Quality controls are carried out for every lot.

Disposal

For any waste produced, each laboratory is responsible for proper disposal in accordance with legal requirements. Be sure to inactivate any infectious material prior to disposal, for example by autoclaving for 30 min. at 134°C.

Reporting information

A report must be submitted to the responsible authority of the Member State in which the user is located and to the manufacturer in the event of a serious incident in connection with this product.

Literature

- ① EUCAST Disk Diffusion Method for Antimicrobial Susceptibility Testing – Version 7.0 (Januar 2019), www.eucast.org
- ② Performance Standards for Antimicrobial Susceptibility Tests; 11th ed., M02-A11; Vol. 32 No.1, Clinical and Laboratory Standards Institute

Abbreviations/Symbols

Symbol	Signification
REF	Catalogue number
IVD	In vitro diagnostic medical device
Σ	Contains sufficient for <n> tests
	Manufacturer
LOT	Batch code
	Do not re-use
CE	CE marking

Change history: Adapted to IVDR specifications



Headquarters

Axon Lab AG
Täferstrasse 15
CH-5405 Baden-Dättwil

www.axonlab.com

Affiliates

Austria
Belgium
Croatia
Czech Republic
Germany

info@axonlab.at
info@axonlab.be
info@axonlab.hr
info@axonlab.cz
info@axonlab.de

Luxembourg
Netherlands
Slovenia
Switzerland (de)
Switzerland (fr)

info@axonlab.lu
info@axonlab.nl
info@axonlab.si
info@axonlab.ch
info_f@axonlab.ch

062022 C 06040007