Utilisation

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Intended use

The iColo-Rectal Test is an immunological screening test for the detection of occult blood in the faeces. This test is used in preventive care of intestinal cancer.

Test principle

The human hemoglobin (hHb) contained in the faeces sample reacts with the monoclonal anti-hHb antibodies (mouse IgG) which are conjugated with colloidal gold particles. Diffusion causes the immune complex to enter the test zone **T** and to react with the immobilized anti-hHb antibodies (goat) to form a purple test line. This denotes a positive test result. The control zone **C** is loaded with immobilized, polyclonal anti-mouse IgG antibodies (goat). By binding with the excess free anti-hHb antibodies (mouse IgG) which are conjugated with colloidal gold particles, these also form a purple line. This is used to validate the proper functioning and implementation of the test.

Reagents/composition

a) The test cassette sealed in aluminum foil (with a desiccant) contains the following:			
Sample well S	Monoclonal anti-hHb antibody gold conjugate (mouse IgG)	15 µg/test	
Test zone T	Polyclonal nitrocellulose-immobilized anti-hHb antibody (goat)	6 µg/test	
Control zone C	Polyclonal nitrocellulose-immobilized anti-mouse IgG antibody (goat)	4 µg/test	
b) Sample tubes PBS buffer, pH 7.5			
		2 ml	

Stability

The reagents are stable up to the indicated expiration date at +4°C to +30°C.

Performing the test

Sample collection


A faeces sample should not be collected or should be collected at a later date given the following:

- During menstruation
- In the case of bleeding hemorrhoids
- In the case of bleeding from constipation because test results would not be correct.

Alcohol, aspirin® and other medications may cause gastrointestinal bleeding and should be discontinued two days prior to faeces sample collection. No specific diet is required prior to performing a test.

Sample collection is performed as follows:

- 1) First collect the stool sample on a clean and dry support (e. g. white paper). The special faeces collecting system (REF 0703380) can also be used to collect stool samples.
- 2) Remove the screw cap of the sample tube and insert the spiral stick approx. 2 cm into the faeces sample to collect a **small quantity of faeces**.
- 3) Insert the spiral stick with the small quantity of faeces into the sample tube with the buffer solution and tighten the screw cap. Shake the tube well so that the faeces loosens from the spiral stick and disperses in the buffer solution.


- 4) Perform this procedure **three times**, i.e., insert the spiral stick into the faeces sample at three different locations and transfer the small quantity of faeces to the buffer solution in the sample tube.
- 5) Write the patient data (name, date of birth, date of sample collection) on the included self-adhesive label and affix the filled-out label in the longitudinal direction to the sample tube.

Convey the sample tube **to the physician or laboratory for analysis as quickly as possible**. The sample may be stored for a maximum of 7 days at +4°C to +30°C or 14 days at +2°C to +8°C.


The procedure for collecting faeces samples is also described in a separate package insert to be handed to the patient with the sample tube incl. label.

Test sequence

- 1) Conduct the test at room temperature (+15°C to +30°C).
- 2) Remove a test cassette from the aluminum foil and place it on a clean, even surface with the sample collector **S** on the right side.
- 3) Shake the sample tube again well and **carefully** break off the tip at the thicker part. Prevent splashing. Place exactly two large drops in the sample well **S** (application point) by squeezing the sample tube.
- 4) A purple line moves slowly to the result window with the letters T and C. After 5 to **max. 10 min.**, the test result is visible as purple lines.

Interpreting the test result

- 1) The purple control line on the left side of the result window (designated by **C**) indicates proper test functioning and implementation.
- 2) The test result is shown by a purple line on the right side of the result window (designated by **T**).

Positive result: The appearance of two color lines (**T** and **C**) denotes a positive result. The intensity of both color lines may vary. A barely visible test line, depending on the hHb concentration, also indicates a positive result. In this case, additional tests by the physician are recommended.

Negative result: The occurrence of only one visible control line **C** indicates a negative result.

Invalid result: If no control lines appear in the result window within 10 min, the result is invalid. The test must be repeated with a new test cassette. An invalid result may occur in case the test components were stored at over 30°C or if not exactly two large drops were placed in the sample well **S**.

This test is a qualitative screening test and therefore cannot be used to definitively diagnose polyps and/or a colorectal carcinoma.

Test characteristics

The detection limit is 50 ng hHb/ml buffer solution (10 µg HB/g feces). The 10 µg HB/g feces cut-off sensitivity is calculated using the method proposed by Fraser et al. to correspond to the cut-off sensitivity of 50 ng/ml buffer solution. In a study involving 260 patients (136 high risk patients), a sensitivity of 82% for colorectal cancer and advanced adenomas > 1cm and a specificity of 99% were determined in comparison to the colonoscopy. The iColo-Rectal Test is specific to human hemoglobin. By adding the following substances to positive and negative samples, no cross reactions or interferences were detected for the indicated concentrations:

Substance	Concentration
Bovine hemoglobin	100 µg/ml
Fish hemoglobin	100 µg/ml
Chicken hemoglobin	100 µg/ml
Horse hemoglobin	100 µg/ml
Sheep hemoglobin	100 µg/ml
Acetylsalicylic acid	20 mg/dl
Ascorbic acid (vitamin C)	20 mg/dl
Glucose	2000 mg/dl
Lithium heparinate	3 mg/dl
Sodium heparinate	3 mg/dl

Hook Effect (Prozone Effect)

The Hook Effect was not observed (false negative results) for Hb concentrations up to 0.1mg/mL.

Precautions and waste disposal

For in vitro diagnostic use.

The analysis material should be considered to be potentially infectious and should therefore be handled with appropriate care.

It is the laboratory's responsibility to comply with the respective national regulatory requirements when disposing of the analysis material and used test components.

References

Fraser_et_al. "A Proposal to Standardize Reporting Units for Fecal Immunochemical Tests for Hemoglobin", Journal of the National Cancer Institute, Oxford Academic, 2012.

Abbreviations/Symbols

Symbol	Signification
	Catalogue number
	In vitro diagnostic medical device
	Contains sufficient for <n> tests
	Manufacturer
	Authorized Representative in the European Community
	Temperature limit
	Use-by date
	Batch code
	Consult instructions for use
	CE Marking
	Do not reuse

® Registered Trade-

	10 Probenröhrchen 10 Testkassetten 10 Etiketten für Probenröhrchen 10 Packungsprospekte für den Patienten
	10 tubes pour l'échantillon 10 cassettes-test 10 étiquettes pour les tubes pour l'échantillon 10 prospectus pour le patient
	10 provette campione 10 cassette per test 10 etichette per provette campione 10 fogli illustrativi per il paziente
	10 sample tubes 10 test cassettes 10 labels for the sample tubes 10 package inserts for the patient

	30 Probenröhrchen 30 Testkassetten 30 Etiketten für Probenröhrchen 30 Packungsprospekte für den Patienten
	30 tubes pour l'échantillon 30 cassettes-test 30 étiquettes pour les tubes pour l'échantillon 30 prospectus pour le patient
	30 provette campione 30 cassette per test 30 etichette per provette campione 30 fogli illustrativi per il paziente
	30 sample tubes 30 test cassettes 30 labels for the sample tubes 30 package inserts for the patient

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© Marque déposée